

Lateral Tongue Localizing (LTL)

- Reversible design allows for left or right tongue localization

Part #	Incisal Separation
0101	1.5 cm
0101B	1.5 cm Slim Paddle


0101B

0101

Mouth Open Tongue Forward (MOTF)

- Upper engagement channel positively engages with either teeth or gums
- Open center section allows the tongue to be placed forward, over the lower incisors
- Lower engagement ridges provide indexing points for the lower jaw

Part #	Incisal Separation
0201	1.5 cm
0202	2 cm
0203	2.5 cm


0201

0203

Mouth Open Tongue Depressed (MOTD)

- Tongue displacement paddle keeps the tongue low, positioned against the base of the mouth
- Upper engagement channel positively engages with either teeth or gums
- Lower engagement ridges provide indexing points for the lower jaw

Part #	Incisal Separation
0300	1.0 cm
0301	1.5 cm
0302	2 cm
0303	2.5 cm


0300

0303

Mouth Open Tongue (Forward) Closed (MOTC)

Part #	Incisal Separation
0401	1.5 cm
0402	2 cm
0403	2.5 cm


0401

0403